CFX-LINKER[®]

corneal cross-linking system

User manual

Rev.11
of 31/03/2023

User manual

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Käyttöopas

CFX-LINKER[®]

corneal cross-linking system

User manual

 English

Rev.11
of 31/03/2023

Table of Contents

Table of Contents 2
Index of figures 3

Section 1 Overview

Overview 4
Information on the user manual 4
Conventions used in this manual 4
Updates to the product and the user manual 4
User comments 5
Information on the device 5
 Control unit 6
 Optical body 6
Information on the device software 7
Labelling 7
Package contents 8
Reproduction 8
Assistance 8

Section 2 Side effects and warnings

Side effects and warnings 9
Side effects 9
Cautions and Precautions 9
Electrical safety 9
Electromagnetic compatibility 10
Patient safety 13
Additional comments on safety 13
Disposal of the product (electrical and electronic waste) 13
Guidelines 14
 Positioning the device 14
 Positioning the operator 14

Section 3 Use of the device

Use of the device 15
Important precautions before switching on the device 15
User Configuration 15
Switching on the device 15
Positioning the device/operator/patient 16
End of treatment 16
Alignment 17

Section 4 Maintenance/Service

Maintenance/Service	18
Installation criteria	18
Maintenance	18
Service	18
Troubleshooting	18
Instructions for sterilisation and disinfection	18
Cleaning the device	19
Cleaning the opening	19
Monitor Cleaning	19
Adjustment of the articulated arm	20
Regular maintenance	20
Handling the device	20
Storage of the device	21
Software	21

Section 5 Features of the device

Features of the device	22
Dimensions of the device	22
Monitor features	22
UV Beam	22
Temperature, humidity and pressure limits	22
Electrical features	22
Reference standards and regulations	23

Index of figures

	Page
Figure 1 Diagram of CF X-LINKER device	5
Figure 2 Diagram of control unit (front)	6
Figure 3 Diagram of control unit (lateral)	6
Figure 4 Optical body diagram	6
Figure 5 Labels affixed to the machine body	7
Figure 6 Labels affixed to both sides of the optical head	7
Figure 7 Diagram showing the positioning of the device/operator/patient	16
Figure 8 System operating diagram for finding the plane of focus	17
Figure 9 Example of correct alignment	17
Figure 10 Example of incorrect alignment	17
Figure 11 Positioning the stand for transportation	20

Section 1 Overview

Read the information in this manual before using the CF X-LINKER medical device.

Information on the user manual

This user manual is a necessary reference tool for use of the device. It is intended for those who already have basic cross-linking skills, and should not be used as a training manual for these skills or other clinical procedures. It is essential to undergo training in cross-linking techniques in order to use the device. The user manual describes the procedures for preparing, using and maintaining the CF X-LINKER medical device. All operating instructions, illustrations, error messages and relevant information are contained in this manual. The operator is responsible for ensuring that all safety instructions found in this manual are applied.

The user manual contains a summary and an analytical index that facilitate the search for information contained therein. The manual is divided into the following sections:

Section 1	“Overview”	General information on the user manual and the device
Section 2	“Side effects and warnings”	Information on safe handling and basic operation practices
Section 3	“Use of the device”	Information on the use of the device
Section 4	“Maintenance/Service”	Information on routine maintenance of the device and the service
Section 5	“Technical features”	Technical features of the device

Conventions used in this manual

The following conventions are used in this user manual.

- The **Caution** notes describe the precautions necessary to avoid causing injury or endangering the life of the user.
- The **Warning** notes describe the precautions necessary to prevent damage to products.
- If the operations in which a procedure is structured are numbered, they should be performed in the corresponding order.
- Bulleted lists present the information in list format, however, this does not imply a sequence.

Updates to the product and the user manual

The Manufacturer will probably offer software updates, new functions and improvements to the device performance in the future. Updates to the user manual accompany software updates and provide detailed information about the improvements made.

User comments

Users are invited to ask questions and make comments on the products by sending these to info@iromedgroup.com.


Information on the device

CF X-LINKER is an electronic medical device used for corneal cross-linking: it emits ultra violet light (wavelength: 365 nm) and enables the photo-polymerisation of stromal collagen fibrils in order to increase the rigidity, resistance and biomechanical stability of the cornea. CF X-Linker is suitable for the most common types of treatment and is equipped with a customisable operation mode, thus enabling “personalised” patient care. The UV emission is controlled by a computer and an advanced software program which, by means of a touch-screen system, enables the selection and execution of the desired treatment.

The CF-X LINKER device is indicated in the conservative treatment of keratoconus and all ectatic corneal diseases. In particular, it is indicated:

- in all ectatic corneal diseases, both genetic and iatrogenic, such as keratoconus and/or marginal pellucid degeneration, a disease characterised by progressive swelling of the cornea and an abnormal laxity of the corneal parenchyma due to reduced cohesion of the collagen lamellae of which it is composed.
- in all ectatic diseases of the cornea caused by infections, trauma, post-LASIK refractive surgery and/or post-photorefractive keratectomy (PRK) and/or radial keratotomy (RK).
- infections in diseases of the cornea, such as keratitis, and generally all situations that require cross-linking.

The following diagram provides a schematic illustration of the device.

Number	Component
1	Control unit
2	Articulated arm
3	Stand
4	Base
5	Optical body

Figure 1 Diagram of CF X-LINKER device

Table 1: Components of the device

Control unit

The control unit is schematically shown below:

Figure 2 Diagram of control unit (front)


Figure 3 Diagram of control unit (lateral)

Number	Component
1	Touch screen monitor
2	Touch screen cable (USB)
3	Power supply
4	Video cable (HDMI)

Table 2: Control unit components

Optical body

The optical body contains the UVA output mechanism. The output area diameter (3-12 mm) can be selected by means of the iris positioned as shown in Figure 4. The purpose of the red lights is to guide the operator during alignment of the device. The operating parameters are entered through a user interface on touch screen.


Figure 4 Optical body diagram

Number	Function
1	Camera connector (Cable 1 optical body)
2	UV system connector
3	UV emission warning light
4	White light
5	LED for identifying the working plane
6	Diaphragm
7	UV Beam

Table 3: Functioning of the optical body

The device works within pre-set values in a range of less than $\pm 20\%$. During operation, the output is constantly monitored (with a frequency greater than one second) by measuring the current supplied to the LED.

Information on the device software

The operation of the CF X-Linker is controlled by the software installed on the device itself. Screen images in this manual are to be considered representative: in fact, there may be no significant differences between what is shown on the device and what is shown in this manual depending on the software version.

Labelling

This paragraph contains the labels located on the instrument and explanations of the symbols used are provided.


Figure 5 Labels affixed to the machine body

Figure 6 Labels affixed to both sides of the optical head

Label	Explanation
	Obligation to consult this manual before use
	Separate collection of waste electrical and electronic equipment (WEEE)
	Warning. Consult instructions for use
SN:	Device serial number
	Year of production
	Manufacturer
CE	CE marking symbol
	Alternating current
	The device could tip over if pushed or not properly attached
	The device emits non-ionising radiation
	Emission of optical radiation
	Relative humidity

Table 4: Explanation of symbols on labelling

Package contents

The dimensions of the CF X-LINKER medical device package are 77 x 57 x 44 cm. Each package contains:

Number	Part
1	Tripod base
2	Wheels
3	Protective cover
4	Optical head in a protective case
5	Control unit
6	Articulated arm
7	Vertical support
8	AC power cord
9	User Manual
10	Software Manual
11	Declaration of Compliance
12	Anti-drying gel
13	Tool set containing screwdriver, USB mass storage device, hex wrench and screws
14	USB stick

Table 5: Package contents

Reproduction

This manual is the exclusive property of the Manufacturer and may not be duplicated, copied or transmitted electronically, even partially, without the written consent from the Manufacturer.

Assistance

In the event of difficulties in using the CF X-LINKER, please contact the authorised area representative.

Section 2 Side effects and warnings

Side effects

This section describes the situations in which the device should not be used, because the risk of use is higher than any possible benefit. Conditions that may contraindicate the use of the device include:

- Cuts or injuries in the eye to be treated
- Chemical damage or delayed epithelial healing in the eye to be treated
- Pregnancy (including during the planning of pregnancy) or breastfeeding
- Corneal thickness less than 300 microns
- Aphakic and pseudophakic patients whose implanted lenses have no UV blocking
- Paediatric patients under 4 years of age.

Cautions and Precautions

- The Custom Fast protocol must only be implemented in combination with the “Ribocross TE” or “Ribo-fast” medical device.
- Assess the potential benefits to patients with the following conditions:
 - Herpes simplex, herpes zoster keratitis, recurrent corneal erosion, corneal dystrophy
 - Epithelial healing disorders

Electrical safety

This equipment requires special precautions regarding electromagnetic compatibility (EMC). It must be installed and used according to the EMC information provided in this manual.

Warning

To avoid the risk of electric shock, this equipment must be connected to a protective earthed socket.

Warning

This equipment operates with hazardous voltages that can cause shocks, burns or lead to death. Do not remove any fixed panels in order to reduce the risk of electric shock, and involuntary exposure to UVA. Make sure that all system services, beyond those described in this manual are only performed by qualified personnel authorized by the Manufacturer.

Caution

Turn off the switch and unplug the power cord before cleaning (disinfection) of the device. Do not pull the cables when removing the power cord from the power outlet. Hold the plug and pull it from the outlet.

Caution

Do not use the device with a damaged power cord.

- Caution**  Place the power cord so that it cannot be stepped on, kinked, pressed, bent, crushed, or accidentally pulled out of the socket.
- Caution**  Do not use the instrument near water, and be careful not to spill liquids on any part of it.
- Caution**  Do not use the CF X-LINKER device in the presence of flammable mixtures or anaesthetics.
- Caution**  Never look directly into the UV light beam. Never direct the light beam on a person except for therapeutic purposes.
- Caution**  Do not use near or stacked with other equipment; if it is used adjacent to or stacked with other equipment, verify that the device operates normally as expected.
- Caution**  The equipment must not be modified.
- Caution**  The use of power cords other than the one supplied is not recommended, since different cables could increase emissions or decrease the immunity of the device.
- Caution**  The mobile terminals may affect the functioning of the device.

Electromagnetic compatibility

The CF X-LINKER device is intended for use in the electromagnetic environment specified below.

Aspects of emission		
Emission test	Compliance	Electromagnetic environment
RF emissions CISPR 11	Group 1	The internal operation of the CF X-LINKER product uses RF energy only. Its RF emissions are therefore very low, and interference with nearby electronic equipment is unlikely.
RF emissions CISPR 11	Class B	The CF X-LINKER product is suitable for use in all buildings, including those used for domestic purposes and those directly connected to a low voltage power supply network that supplies domestic buildings.
Harmonic emissions IEC 61000-3-2	Class A	The appliance can be used in all buildings, including domestic buildings, and those directly connected to the public low voltage power supply network that supplies buildings for domestic use.
Emission of voltage fluctuations/flickers IEC 61000-3-3	Compliant	

Aspects of immunity

The CF X-LINKER product is intended for use in the electromagnetic environment specified below. The client or user should ensure that it is used in such an environment

Immunity tests	IEC 60601 test level	Level of compliance	Electromagnetic environment
Electrostatic discharge (ESD) IEC 61000-4-2	± 8kV on contact ±15kV in air	± 8kV on contact ±15kV in air	Floors should be made of wood, concrete or ceramic tile. If the floors are covered with synthetic materials, the relative humidity should be at least 30%.
Transients/rapid electrical pulse sequence IEC 61000-4-4	± 2kV for power supply lines	± 2kV for power supply lines	The quality of the mains voltage should be that of a typical commercial or hospital environment.
Over voltages IEC 61000-4-5	±1/2kV differential mode	±1/2kV differential mode	The quality of the mains voltage should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage fluctuations on power input lines IEC 61000-4-11	<5% UT (>95% gap in UT) for 0.5 cycles <5% UT (>95% gap in UT) for 1 cycle 70% UT (30% gap in UT) for 25 cycles <5% UT (>95% gap in UT) for 250 cycles	<5% UT (>95% gap in UT) for 0.5 cycles <5% UT (>95% gap in UT) for 1 cycle 70% UT (30% gap in UT) for 25 cycles <5% UT (>95% gap in UT) for 250 cycles	The quality of the mains voltage should be that of a typical commercial or hospital environment. If the user requires continued operation during mains voltage interruptions, it is recommended that the appliance is powered by an uninterruptible power supply or batteries.
Magnetic field at mains frequency (50Hz) IEC 61000-4-8	30A/m	30A/m	Magnetic fields at network frequency should be at levels that are normal for a typical location in a commercial or hospital setting.

NOTE: UT is the AC mains voltage before applying the standard of proof.

Aspects of immunity to R.f.

The CF X-LINKER product is intended for use in the electromagnetic environment specified below. The customer or user should ensure that it is used in such an environment.

Immunity test	Standard of proof EN 60601-1-2	Level of compliance	Electromagnetic environment GUIDE
RF conducted EN 61000-4-6	6Veff from 150kHz to 80MHz	6Veff from 150kHz to 80MHz	Portable and mobile Rf communication devices should not be used near any part of the device, including cables, unless they comply with the recommended separation distances calculated using the applicable equation for the transmitter frequency. Recommended separation distances: d = 1.2 x √P from 150kHz to 80MHz d = 1.2 x √P from 80MHz to 800MHz d = 2.3 x √P from 800MHz to 2.7GHz Where P is the maximum output power of the transmitter in Watts (W) according to the transmitter manufacturer, and (d) is the recommended separation distance in metres (m)
RF Radiated EN 61000-4-3	10Veff from 80MHz to 2.7GHz	10Veff from 80MHz to 2.7GHz	

The magnetic field strength of fixed RF transmitters, as determined in an electromagnetic site survey, may be less than the compliance level in each frequency range.

Interference may occur near devices marked with the special symbol 

The device is intended for use in an electromagnetic environment in which radiated radio frequency disturbances are controlled. The patient or family member using the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, in connection with the peak output of the communications equipment.

Maximum output power of the specified transmitter (W)	Separation distance to the transmitter frequency (m)			
	From 15 kHz to 80MHz	From 80MHz to 800MHz	From 800MHz to 1GHz	From 1GHz to 2.5GHz
0.01	0.12	0.04	0.07	0.23
0.1	0.37	0.11	0.22	0.74
1	1.17	0.35	0.70	2.33
10	3.69	1.11	2.21	7.38
100	11.67	3.50	7.00	23.33

Patient safety

The treatment must take place in a calm and relaxed environment so as not to distract the patient's attention. The patient must lie on a couch. The patient's head should rest comfortably on a headrest.

Caution

It is imperative that neither the couch nor the CF X-LINKER device are moved during treatment.

Caution

The CF X-LINKER is a medical device. It may only be used in healthcare facilities or medical areas, by medical specialists and practitioners (e.g. optometrists), within the limits permitted by law and according to regulations concerning professional practice.

Additional comments on safety

Caution

The UV beam is calibrated in the factory and is checked once a year. **Any external modifications** (for example, interposition of optical elements) are prohibited.

Caution

Any plastic equipment can be damaged if affected by UV beam resulting in possible product degradation.

Caution

It is recommended to wear protective goggles when using the device to avoid potentially hazardous light radiation exposure.

Caution

Before use, place the device so that it is easier to disconnect in case of need.

Caution

Do not use in oxygen-rich environments.

Caution

The device is classified under risk class 2 (EN 60601-2-57: 2011) up to a distance of 200 cm from the emitter of the UV beam.

Caution

Avoid unauthorised use of the device.

Disposal of the product (electrical and electronic waste)

This marking shown on the product indicates that the product must not be disposed of with other household waste at the end of the life cycle. To prevent possible damage to the environment or health from uncontrolled waste disposal, the user is encouraged to separate the product from other types of waste and recycle them responsibly to promote the sustainable reuse of material resources. Users should contact their supplier and check the terms and conditions of the purchase contract. Do not mix this product or its electronic accessories with other commercial waste.

Guidelines

The purpose of the guidelines set out below is to enable the practitioner to use the device as comfortably and efficiently as possible.

Caution

Use of the device refers to the physical interaction between the practitioner and the device. As with other physical activities of a similar nature, the user may occasionally experience discomfort in the hands, fingers, arms, shoulders, eyes, back or other parts of the body when using the CF X-LINKER device. However, symptoms such as constant or recurring discomfort, pain, spasms, intense pain, tingling, numbness, burning sensations, or stiffness should not be underestimated and medical advice should be sought. These symptoms can be linked with musculoskeletal disorders (MSD), which, besides being painful, can result in potentially disabling injuries to the nerves, muscles, tendons, or other body parts. Examples of MSD include tunnel carpal tunnel and tendinitis. Researchers are not able to give definitive answers to the many questions related to MSD, however, it is generally agreed that the factors associated with their occurrence include: pre-existing medical and physical conditions, overall health conditions, positioning of the equipment and body when working, working frequency, hours of work and other at-risk physical activities. This chapter provides guidelines that may help improve comfort and reduce the risk of developing MSDs..

Positioning the device

To reduce eye strain to a minimum, take the following precautions:

- Position the device within the visual radius.
- Adjust the angle of the device or display to minimise the reflection of artificial or external lighting.

To reduce neck strain to a minimum, take the following precautions:

- Position yourself so that the screen is slightly under eye level.

Positioning the operator

To support your back, take the following precautions:

- Use a chair with lumbar support
- The chair must be adjusted based on the height of the work surface and foster a natural body posture.
- The chair must enable rapid height adjustment.
- Always sit or stand with your back straight. Avoid bending or curving the back.

To reduce extensions and twists to a minimum, take the following precautions:

- Use a couch with an adjustable height.
- Get as close to the patient as you can.
- Face forward. Avoid turning the head or body.

Schedule breaks:

- Minimising working time and scheduling breaks can be very useful to compensate for the associated physical activity and avoid stress, therefore, the onset of MSD.
- Stopping and relaxing may be enough, however, simply changing tasks can help some muscle groups relax while others remain or become active.

Section 3 Use of the device

1. Position the CF X-LINKER device close to the couch. Lock the wheels to ensure that the device is stable.
2. Make sure that the two cables are correctly connected to the optical body according to the colour coding.
3. Make sure that the power cable of the optical body is connected.
4. Make sure there is no dirt or dust on the UV beam opening window.

Important precautions before switching on the device

Caution

The operator must ensure that the CF X-LINKER device is functioning correctly and is in good condition before the treatment begins.

User Configuration

The X-Linker CF provides three types of user accounts: **administrator**, **operator** and **technical service**. If the privacy configuration is activated, the user login will be required to access the Archive section, the Control Panel (administrator and technical service only) and the section on filling in the Treatment Report. Only the administrator can operate on the user database, adding and/or deleting accounts (please refer to the System Administrator's Guide). Only the technical service account can access the Control Panel to proceed with Installation and Maintenance operations.

Switching on the device

1. Connect the AC cable to the power socket.
2. To switch on the device, press the red switch located on the base.
3. Then press the "ON" switch. This will turn on all system components.
4. The CF X-LINKER device will launch the sequence for powering up, loading the operating system and self-configuration.

Caution

If an error occurs at this stage, make a note of the messages issued by the device and contact the distributor or the technical assistance service immediately.

For correct use of the device please refer to the User Manual – Software.

Positioning the device/operator/patient

To execute treatment, we recommend the positioning showed in image:


Figure 7 Diagram showing the positioning of the device/operator/patient

Specifically:

- Make sure the patient is lying flat on a table or on a cot, and that his/her head is steady.
- Adjust the couch so that the patient can rest comfortably for the entire duration of the treatment without moving his/her head.
- Attach any accessories required for the treatment and/or required by the clinical standards.
- Apply the riboflavin-based solution on the area to be treated, according to the protocol.

Caution


Make sure the patient and the CF X-LINKER device are kept still and do not move after alignment or during treatment.

End of treatment


The software alerts the practitioner when the treatment is complete, enabling them to carry out further treatment or switch off the device.

Alignment


The device is equipped with a system that enables the operator to identify the working surface. When required, the device emits two beams of red light that converge at a point located on the plane of focus.

Figure 8 System operating diagram for finding the plane of focus

The plane of focus is correctly detected if the operator sees the following image:

Figure 9 Example of correct alignment

The plane of focus **is not correctly detected** if the operator sees an image similar to the following:

Figure 10 Example of incorrect alignment

Section 4 Maintenance/Service

In this manual, “maintenance” refers to a series of non-technical procedures that must be carried out by the operator in order to keep the device in optimal working order. On the other hand, “service” refers to a series of activities that must be carried out exclusively by a qualified technician.

Installation criteria

For each new client, the authorised distributor provides a complete initial installation and start-up of the device. Following initial installation, and once the device is functioning correctly, the representative may also offer basic training for a designated practitioner in the operation of the CF X-LINKER device. For this reason, no specific instructions for installing or setting up the device are included in this manual. Under the service contract, any additional hardware adjustment, other than what is specified for normal operations, must be performed by (or under the guidance of) a distributor authorised by the Manufacturer.

Maintenance

The client does not generally need to carry out maintenance on the CF X-LINKER device. All maintenance or technical assistance interventions will be carried out by a qualified technician. In the event of problems with the device, please contact the local representative designated by the Manufacturer.

Service

A Service contract is provided with all CF X-LINKER devices; this includes scheduled maintenance and any updates. Moreover, the possibility of unplanned calls is included if required. Please refer to the Manufacturer’s designated local representative.

Troubleshooting

The CF X-LINKER device includes an initial calibration check to confirm that the device is functioning correctly. If the check is not successful, the software will not be able to start a treatment. In this case, the check may be repeated. If this is also unsuccessful, authorised assistance must be contacted.

Instructions for sterilisation and disinfection

None of the components of the CF X-LINKER device require sterilisation: only external cleaning and disinfection are recommended. Only use isopropyl alcohol sprays or preparations for disinfection. Use small amounts of liquids and soft wipes without fibre.

Cleaning the device

- Use a soft, damp cloth.
- The outside of the CF X-LINKER device can be cleaned using a lint-free cloth moistened with diluted bleach, soap and water, or isopropyl alcohol.
- If necessary, a 70% alcohol or 10% bleach solution can be used.
- Do not immerse the device in liquids or pour liquids over the device.
- The glass of the beam opening must not enter into contact with any of the aforementioned solutions for whatever purpose.
- Make sure that liquid detergents do not leak into the device while the outer surfaces are being cleaned. Any leaks could damage the device.

Caution

Disconnect the power supply before cleaning the device.

Caution

Harsh detergents, especially those containing abrasives or corrosive solvents can damage the surfaces of the components.

Cleaning the opening

- Check that the UV beam opening is clean before treating a patient.
- Use wipes specific to camera lenses or compressed air to remove dust and particles from the opening glass surface.

Monitor Cleaning

- Never clean the screen when it is switched on or connected. Switch off the monitor and disconnect the power cord and input cables.
- Disconnect all devices connected to the monitor.
- Do not use liquids directly on the screen.
- Start with a slightly wet microfibre cloth. The cloth should be wet, but not completely soaked.
- If water alone does not work, a specialised screen cleaning product may be used. However, such a product should not contain benzene, solvents, ammonia, abrasive materials, alcohol or acetone. Do not allow chemical agents to come into direct contact with the monitor. Such substances may visibly damage the screen.
- Do not apply too much pressure. High pressure can dull the pixels, causing image uniformity problems, and also brightness problems. Do not use brushes, as they may cause abrasions and scratches in the panel.
- Ensure that surfaces are completely dry before switching on the monitor. No liquid or moisture should be visible in the surface of the product.

Adjustment of the articulated arm

If the arm does not keep the optical head in a fixed vertical position, immediately contact a dealer authorised by the Manufacturer.

Regular maintenance

The calibration status must be checked at least once a year by specialist staff authorised by the Manufacturer. If necessary, they will re-calibrate the device.

Handling the device

The CF X-LINKER device can be moved around within the healthcare facility. The Manufacturer is not responsible for any damage to the device caused by transferring it from one building to another. Before moving the device from one room to another, make sure that the articulated arm is positioned as shown in figure:

Figure 11 Positioning the stand for transportation

Once placed in this position, the device can easily be transported by pushing the stand.

Storage of the device

- Follow the temperature and humidity specifications stated in the respective section.
- Switch off all components and disconnect the power supply.
- Cover the touch screen display on the control unit.
- To avoid damage or misalignment, do not tamper with any parts of the device.

Software

If there are errors or malfunctioning of the software, immediately contact the Manufacturer or one of the area representatives.

Caution

Visually inspect the device before use. Do not use it if damage or malfunctioning are evident. Using a device that is not functioning correctly may cause harm to the operator and/or patient.

Section 5 Features of the device

CF X-LINKER is a UVA illuminator to be used on the cornea: it provides light radiation of varying intensity and duration at a nominal fixed wavelength of 365 nm. It is used for therapeutic and/or preventive purposes. As stipulated in Annex IX of Directive 93/42/EEC and subsequent amendments. CF X-LINKER is a class IIa medical device according to rule 9.

Dimensions of the device

Base diameter: 60 cm

Height: 126 cm

Maximum arm extension: 90 cm

Weight: Approximately 30 kg

Monitor features

Diagonal: 254 mm (10 inches)

Resolution: 1920 x 1200 pixels

UV Beam

Wavelength: 365 nm

Intensity: 1 to 30 mW/cm²

Mode: Continuous and Pulsed

Diaphragm opening: 3-12 mm

Protocols: Dresden, Custom Fast, Iontophoresis, Keratitis, Preset, Demo

Temperature, humidity and pressure limits

The device must be used at a temperature of +0 to +40 °C, transported and stored at a temperature of 0 to +70 °C. It must also be stored, transported and used in a relative humidity of less than 70% and an atmospheric pressure of 500 to 1060 mbar.

Electrical features

Input power: 100 - 240 V AC.

Reference standards and regulations

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Iromed Group Srl