



Innovation

Rx  
The first  
**FEMTOLASER  
INSTRUMENT**  
for capsulotomy



*Give your  
advanced IOLs  
the best*

# Rx, the first FEMTOLASER INSTRUMENT for capsulotomy



Non-contractual pictures

## Mobile and ultra compact

Predictable and safe

Battery powered

Handheld instrument

Mobile

Compact (50x50x110 cm)




## Handheld instrument: an ergonomic handpiece




CE marking pending

## Round, centered and predictable capsulotomy



Rx performs the anterior capsulotomy by a Helix-shaped laser beam cut:

- Capsulotomy: Ø 5 mm,
- Procedure: 4 s,
- Bag cut: 20 ms



ILASIS LASER  
8 rue Adrienne Bolland, 33600 Pessac  
+33 (0)5 64 31 19 30 - contact@ilasis.com



[www.helixsurgical.com](http://www.helixsurgical.com)

Scan to watch



RX\_FLYER\_ED1\_EN\_11\_2023