



The future of vision

SIOBAL








SIOBAL Silicone Oil



<i>Brand name</i>	SIOBAL		
<i>Reference</i>	S1000	S2000	S5000
<i>Description and indications</i>	High-purity silicone oil, free of low-molecular-weight components. Tamponade agent for retinal detachment. <ul style="list-style-type: none"> • Retinal detachment • Proliferative vitreoretinopathy. • Perforating wounds. • Proliferative diabetic retinopathy. 		
<i>Viscosity (1s⁻¹; 25°C)</i>	1000 cSt	2 000 cSt	5 000 cSt
<i>Sterilization method</i>	Steam.		
<i>Supply</i>	Sterile, single unit. 10ml prefilled glass syringe.		
<i>Shelf-life</i>	3 years		





* Product with CE marking

SIOBAL PACK

Brand name	SIOBAL PACK: SIOBAL Silicone Oil + Injection Tubing Set (Sterile)											
References	PS1000 *				PS2000 *				PS5000 *			
Models	A	OA	OR	B	A	OA	OR	B	A	OA	OR	B
	*Selection of compatible connection system for the following equipments: (A) Alcon: Accurus, Constellation Optikon: (OA) Assistant, (OR) R-Evolution. (B) Bausch&Lomb: Millenium, Stellaris											
Description and indications	SIOBAL silicone oil injection /extraction system. This system is connected to vitrectomy equipment so as to facilitate the SIOBAL silicone oil injection/extraction process.											
Viscosity- oil	1000 cSt				2000 cSt				5000 cSt			
Content	Cannula 20G	Cannula 23G	Cannula 25G	Connection system for injection / extraction				10ml oil prefilled syringe				
												
Supply	5 packs per box. Products sterilized and packed in their original containers by the manufacturer.											
Quality assurance	Products with CE marking											



SIOBAL SET

Brand name	SIOBAL SET: Silicone Oil Injection / Extraction Tubing Set (Sterile)				
Reference	ST *				
Models	A	OA	OR	B	
	*Selection of compatible connection system for the following equipments: (A) Alcon: Accurus, Constellation Optikon: (OA) Assistant, (OR) R-Evolution. (B) Bausch&Lomb: Millenium, Stellaris				
Content	Cannula 20G	Cannula 23G	Cannula 25G	Connection system for injection / extraction	
					
Supply	5 sets per box. Products sterilized and packed in their original containers by the manufacturer.				
Quality assurance	Products with CE marking				

REV: 04. 06/2019