


Sulco Fix Add-on

Minimum Aberration, Maximum Acuity



Available in
Monofocal, Toric, Trifocal and Toric Trifocal

Pseudophakic Add-ON IOLs



- Excellent Quality of Vision.
- Fast Recovery.
- Very Stable Over the time.
- Customized Optic and Haptic Profile

Available in
Monofocal, Toric, Trifocal and Toric Trifocal

Specification:	Model	UAON	UAONT	UAOND	UAONDT
Overall Length	Material	Hydrophilic Acrylic			
	IOL Powers (Spherical Equivalent Diopters)	-15.00 to +15.00 D (in 0.5 D Steps)			
	IOL Cylinder Power	+0.50 D to +10.00 D (in 0.5 D Steps)			
	Addition	+3.5 D for Near +1.8 D for Intermediate			
	Optic Diameter	6.20 mm			
	Overall Length	11.00 to 14.00 mm (in 0.25 mm Steps)			
	Optic Type	Aspheric - Convex/Concave			
	Vault Height	0.00 mm			
	Refractive Index	1.465			



CAREGROUP SIGHT SOLUTION LLP
 Block No.: 310/C&E, Village Dabhasa, Taluka - Padra,
 Dist. Vadodara - 391 440, Gujarat, India.
 Customercare@caregroupindia.com
 www.caregroupindia.com | www.caregroupiol.com