



Basis Z

Technische Spezifikationen



hydrophob preloaded 


A-Konstante* 119,033 (SRK/T)	Holladay 1* 1,809 (SF)
Haigis Konstanten* a0 = 0,8751 a1 = 0,2659 a2 = 0,1416	Holladay 2 5,57 (ACD)
Hoffer Q* 5,60 (pACD)	Barrett* 1,90 (LF); N/A (DF)
	*Stand IOL Con 03.2023

Refraktiver Index 1,47	Abbe-Zahl 58
----------------------------------	------------------------

Linsentyp 1-teilige Hinterkammerlinse	Sterilisation EO
Material biokompatibles hydrophobes Acryl (<0,5%)	Empf. Injektionssystem POB-MA
Haptik Z-Haptik, 0° Anwinkelung	Filter UV-Filter, Blaulichtfilter (optional)
Optik asphärische Optik, 360° Zellbarriere	

Lieferbereich

Artikelnummer klar	Artikelnummer gelb	Sphäre	Schritte
S1AB00	S1ABY0	-10,0 D bis -1,0 D	1,0 D
B1AB00	B1ABY0	0,0 D bis +9,0 D	1,0 D
B1AB00	B1ABY0	+10,0 D bis +30,0 D	0,5 D
S1AB00	S1ABY0	+31,0 D bis +35,0 D	1,0 D



Immer die aktuelle Version. QR-Code scannen.
www.1stq.de/basisz1021



Haben Sie noch Fragen?
+49 621 7176330 • info@1stq.de

